

SIUI

SHANTOU INSTITUTE OF ULTRASONIC INSTRUMENTS

CTS-3300V

Digital Veterinary Ultrasound Imaging System

- ④ 10-inch High Resolution Monitor
- ④ Broadband Multi-frequency Probes
- ④ Non-volatile Storage Capacity
- ④ Veterinary Software Package



CTS-3300V is a new portable ultrasound imaging system, which can be widely applied to various vet diagnosis. It has non-volatile storage capacity and abundant veterinary measurement and calculation package. With superb image quality and stable performance, it surely will bring you a fresh and new experience in ultrasound exam.

Standard Configuration:

- Main Unit 1 Unit
- Linear Probe 1 Pc
- Accessories 1 Set

Options:

- 2.5~5.0MHz Convex Probe
- 3.5~7.5MHz Vaginal Probe
- B/W Video Printer
- Biopsy Guide
- Foot switch
- Trolley

Multi-frequency




Image transmission



Images can be transferred to PC through RS-232C interface to save and process in JPG, BMP or TIF format, enabling exchange of medical images between medical devices on a network, easy to store and print out.

Veterinary OB Software:

BOVINE

- BPD (Biparietal Diameter)
- BD (Body Diameter)
- CRL (Crown-Rump Length)
- TD (Trunk Diameter)
- VOL-B (Volume)
- ANGLE (Angle)

CANINE

- GSD (Gestational Sac Diameter)
- CRL (Crown-Rump Length)
- HD (Head Diameter)
- BD (Body Diameter)
- VOL-B (Volume)
- ANGLE (Angle)

FELINE

- HD (Head Diameter)
- BD (Body Diameter)
- VOL-B (Volume)


OVINE

- BPD (Biparietal Diameter)
- USD (Umbilicus to Spine Diameter)
- CRL (Crown-Rump Length)
- TD (Trunk Diameter)
- VOL-B (Volume)
- ANGLE (Angle)

EQUINE

- BPD (Biparietal Diameter)
- GSD (Gestational Sac Diameter)
- VOL-B (Volume)
- ANGLE (Angle)

High-quality



Features:



Backlit Keyboard



10-inch High Resolution Monitor



Trolley

SHANTOU INSTITUTE OF ULTRASONIC INSTRUMENTS

Add:#77, Jinsha Road, Shantou 515041, Guangdong, China
 Tel: (86)754-8250150 Fax: (86)754-8251499
 Http://www.siui.com E-mail:siui@siui.com

