

# Midmark 646 & 647

Power Podiatry Procedures Chair  
Foot, Ankle & Wound Care Procedures



**Midmark 647-002**  
shown in Silver Sage  
Premium seamless  
upholstery

The **Midmark 646 and 647 Barrier-Free® Power Podiatry Procedures Chairs** are the next step in foot, ankle and wound care. The new, ergonomic foot section design gives you what you need at the point of care to perform even the lengthiest procedures, with many standard features:

- Convenient, integrated touch controls on both sides of the chair
- Seamless, easy-to-clean surface
- Barrier-Free® low height of 19 inches improves safety
- Fast, efficient movement
- Programmable and rotation models available
- Wide patient support surface of 28 inches accomodates all patients

**Midmark 647-002**  
shown in Ultra-Arctic  
Ultra-Premium plush  
upholstery (additional  
cost: \$289)



## Midmark 646 & 647 Features



### Accessibility

A barrier-free low height of 19" and standard, adjustable chair arms make transfer a breeze for all patients.



### Release Handles

Ergonomic, conveniently positioned release handles allow for smooth, one-handed operation (patent pending).



### Adjustable Arms

Standard, adjustable chair arms rotate completely out of the way for unsurpassed comfort and accessibility for all patients.



### Debris Tray

Large stainless steel debris tray with no exposed glides.

# Promotional Offers

For each Midmark 646 or 647 purchased between **March 1 – May 31, 2010**, with delivery taken by **August 31, 2010**, you will receive your choice of **one** of the following from Midmark:

Vision Block Screen

- AND -

Ritter 276 Air Lift Stool

Total package valued at

**\$818**

**\$400**

Rebate from Midmark

Four (4) year warranty from Midmark

(three (3) year extended + one (1) year standard warranty)

Valued at

**\$1,200**

## Eligible Products & Pricing

	List Price	Promo Price	Savings
<b>Midmark 646-001</b> (non-programmable with receptacles, fixed base, premium upholstery)	\$6,212	\$5,079	<b>\$1,133</b>
<b>Midmark 647-001</b> (non-programmable with receptacles; premium upholstery)	\$8,493	\$6,807	<b>\$1,686</b>
<b>Midmark 647-002</b> (programmable with receptacles; premium upholstery)	\$9,036	\$7,507	<b>\$1,529</b>
<b>Midmark 647-003</b> (non-programmable with rotation, without receptacles; premium upholstery)	\$10,102	\$8,407	<b>\$1,695</b>
<b>Midmark 647-004</b> (programmable with rotation, without receptacles; premium upholstery)	\$10,783	\$9,107	<b>\$1,676</b>



### Covered Glides

Covered glides to simplify cleaning and maintenance (patent pending).



### Touch Controls

Integrated touch controls on both sides eliminates costly hand control breakage and allows unrestricted patient access.



### Rotation (model specific)

Improve efficiency by positioning your patient instead of peripheral equipment. The simple touch of a pedal offers 350 degrees of rotation.

# Product Comparison Chart

Features	641-002	641-003	641-004	641-005	630-003	630-004	630-005	630-006	646-001	647-001	647-002	647-003	647-004	230-001	230-002	230-003	244-001
Patient weight capacity	450 lbs.	450 lbs.	450 lbs.	450 lbs.	450 lbs.	450 lbs.	450 lbs.	450 lbs.	450 lbs.	450 lbs.	450 lbs.	450 lbs.	450 lbs.	450 lbs.	450 lbs.	450 lbs.	850 lbs.
Minimum seat height	19"	19"	22.5"	22.5"	19"	19"	22.5"	22.5"	21"	19"	19"	22.5"	22.5"	22.5"	22.5"	26"	18"
Maximum seat height	40"	40"	43.5"	43.5"	40"	40"	43.5"	43.5"	21"	40"	40"	43.5"	43.5"	40"	40"	43.5"	34"
Headrest upholstery dimension	9"	9"	9"	9"	12"	12"	12"	12"						18"	18"	18"	
Back upholstery dimension	21.5"	21.5"	21.5"	21.5"	24"- 26"	24"- 26"	24"- 26"	24"- 26"	28"	28"	28"	28"	28"	27"	27"	27"	32"
Seat upholstery dimension	23.5"	23.5"	23.5"	23.5"	28"	28"	28"	28"	28"	28"	28"	28"	28"	27"	27"	27"	32"
Foot upholstery dimension	21.75"	21.75"	21.75"	21.75"	18"	18"	18"	18"	26"	26"	26"	26"	26"	18"	18"	18"	30"
Foot control	•	•	•	•	•	•	•	•	•	•	•	•	•	• opt	• opt	• opt	• opt
Hand control	•	•	•	•	•	•	•	•	← touch controls →					•	•	•	•
QuickChair™/Auto-return	•	•	•	•	•	•	•	•	•	•	•	•	•				
Crash-avoidance safety system	•	•	•	•	•	•	•	•		•	•	•	•				
Premium Comfort System™	•	•	•	•	•	•	•	•	•	•	•	•	•				
Hospital grade electrical outlets		•			•	•			•	•	•				•		
Side rail adjustment	•	•	•	•													
Flip-out foot extension	•	•	•	•					← slide out extension →								
Stirrups/Stainless steel treatment pan					•	•	•	•	← stainless steel pan →					•	•	•	
Knee-chest positioning														•	•	•	
Programmability		•		•		•		•			•		•				
Rotation			•	•			•	•				•	•			•	
Paper roll holder	• opt	• opt	• opt	• opt	•	•	•	•		• opt	• opt	• opt	• opt	•	•	•	•
Wireless controls	• opt	• opt	• opt	• opt	• opt	• opt	• opt	• opt		• opt	• opt	• opt	• opt				
Ultra-Premium upholstery	• opt	• opt	• opt	• opt	• opt	• opt	• opt	• opt		• opt	• opt	• opt	• opt				
Upholstery colors	27	27	27	27	27	27	27	27	27	27	27	27	27	8	8	8	8

opt = optional

## Midmark Premium Solid Upholstery Colors



Blueberry -235



Cashmere Blue -239



Silver Sage -236



Terra Cotta -240



Fossil Grey -238



Tea Green -237



Black -312

## Midmark Premium Pattern Upholstery Colors



Belagio -247



Capri -251



Siena -248



Treviso -252



Firenze -250



Tuscany -249

## Midmark Ultra-Premium Solid Upholstery Colors



Ultra-Pacific Blue -241



Ultra-Sky -245



Ultra-Spruce -242



Ultra-Canyon -246



Ultra-Arctic -244



Ultra-Olivine -243



Ultra-Raven Wing -228

## Ritter Solid Upholstery Colors



Navy -231



Dusty Blue -233



Moss -230



Clay -234



Pebble Grey -216



Shadow -232



Perfect Plum -229



Black -312

For more information or a demonstration, contact your Midmark or Ritter Dealer. Midmark Corporation, 60 Vista Drive, P.O. Box 286, Versailles, Ohio 45380-0286. Call: 1-800-MIDMARK Fax: 1-800-365-8631 midmark.com